# Rapid SARS-CoV-2 Antigen Test Card Self-Test

CE  
0123

**Xiamen Boson Biotech Co., Ltd.**

**Product Brochure**

## Rapid SARS-CoV-2 Antigen Test Card

- Self-test received **notified body approved CE marking** on April 1, 2021
- A **home test** for the rapid qualitative detection of SARS-CoV-2 virus antigen in anterior nasal swabs within 7 days of symptom onset

# Product Information

## Performance Characteristics

- **Sample type:** nasal swab
- **Time to result:** 15-20 min
- **Sensitivity:** 96.77%
- **Specificity:** 99.20%
- **Accuracy:** 98.72%
- **Layperson study statistics:**
  - 99.10% of non-professionals carried out the test without requiring assistance
  - 97.87% of different types of results were interpreted correctly

## Specifications

- 1 Test / Kit  
 Catalog No.: 1N40C5-2  
 Kit Components:
- 1 Rapid SARS-CoV-2 Antigen Test Card (sealed foil pouch)
  - 1 sterilized swab
  - 1 extraction tube
  - 1 sample extraction buffer
  - IFU
  - Tube holder on back of box





- 5 Tests / Kit  
 Catalog No.: 1N40C5-4  
 Kit Components:
- 5 Rapid SARS-CoV-2 Antigen Test Card (sealed foil pouch)
  - 5 sterilized swab
  - 5 extraction tube
  - 5 sample extraction buffer
  - IFU
  - 1 tube stand

- 20 Tests / Kit  
 Catalog No.: 1N40C5-6  
 Kit Components:
- 20 Rapid SARS-CoV-2 Antigen Test Card (sealed foil pouch)
  - 20 sterilized swab
  - 20 extraction tube
  - 20 sample extraction buffer
  - IFU
  - 1 tube stand

## Test Procedures

①		<p>Rotate the lid of sample extraction buffer bottle. <b>Caution: Open it away from your face and be careful not to spill any of the liquid.</b></p>
②		<p>Squeeze all extraction buffer out of the bottle into the extraction tube. <b>Caution: Avoid touching the bottle against the tube.</b></p>
③		<p>Identify the soft, fabric tip of the swab. Peel open the swab packaging and gently take out the swab. <b>Caution: Never touch the soft, fabric tip of the swab with your hands.</b></p>
④		<p>Carefully insert swab into one nostril. The swab tip should be inserted no less than 2.5 cm (1 inch) from the edge of the nostril. Roll swab 3-4 times along the mucosa inside the nostril. Leave swab in the nostril for several seconds. Using the same swab, repeat this process for the other nostril. Withdraw swab from the nasal cavity. <b>Caution: This may feel uncomfortable. Do not insert the swab any deeper if you feel strong resistance or pain.</b></p>
⑤		<p>Place swab into extraction tube. Roll swab three to five (3-5) times. <b>Leave swab in extraction buffer for 1 minute.</b></p>
⑥		<p>Pinch extraction tube with fingers and remove the solution from swab as much as possible.</p>

⑦		Install the nozzle cap onto the sample extraction tube tightly.
⑧		Bring the kit components to room temperature before testing. Open the pouch and remove the card. Place the card on a flat and level surface. <b>Caution: Once opened, the test card must be used immediately.</b>
⑨		Invert the extraction tube and add 3 drops of test specimen into the specimen well (S) by gently squeezing the extraction tube. <b>Caution: The formation of air bubbles in the specimen well (S) must be avoided.</b>
⑩		Read the results at 15-20 minutes. <b>Caution: Results after 20 minutes may not be accurate.</b>

## Interpretation of Results


	<p><b>Positive:</b> If two colored bands appear with one colored band in the Control Zone (C) and another in the Test Zone (T) within 15-20 minutes, the test result is positive.</p> <p><b>Caution:</b> No matter how faint the colored band is in the Test Zone (T), the result should be considered as positive.</p>
	<p><b>Negative:</b> If one colored band appears in the Control Zone (C) and no colored band appears in the Test Zone (T) within 15-20 minutes, the test result is negative.</p>
	<p><b>Invalid:</b> If no color line appears in the control area (C) within 15-20 minutes, the test is invalid. Repeat the test with a new test card.</p>

# Authorizations and Certificates

# CE Certificate by Notified Body

ZERTIFIKAT ◆ CERTIFICATE CERTIFICADO ◆ CERTIFICAT CERTIFICATE ◆ 認證證書		Benannt durch/Designated by Zentralstelle der Länder für Gesundheitsschutz bei Arzneimitteln und Medizinprodukten www.zlg.de ZLG-BS-245.10.07		
	Product Service			
	<b>EC Certificate</b>			
	EC Design-Examination Certificate Directive 98/79/EC on In Vitro Diagnostic Medical Devices (IVDD), Annex III (6) (Devices for self-testing)			
	<b>No. V9 061317 0006 Rev. 00</b>			
	<b>Manufacturer:</b>		<b>Xiamen Boson Biotech Co., Ltd.</b> 90-94 Tianfeng Road Jimei North Industrial Park 361021 Xiamen, Fujian PEOPLE'S REPUBLIC OF CHINA	
	<b>Product:</b>		In Vitro diagnostic devices for self testing	
	<b>Model(s):</b>		Rapid SARS-CoV-2 Antigen Test Card	
	<b>Parameters:</b>		Model Name:	Model No.:
			Rapid SARS-CoV-2 Antigen Test Card	REF 1N40C5-2 REF 1N40C5-4 REF 1N40C5-6
The Certification Body of TÜV SÜD Product Service GmbH declares that a design examination has been carried out on the respective devices in accordance with IVDD Annex III (6). The design of the devices conforms to the requirements of this Directive. All applicable requirements of the testing and certification regulation of TÜV SÜD Group have to be complied with. For details and certificate validity see: <a href="http://www.tuvsud.com/ps-cert?q=cert.V9.061317.0006.Rev.00">www.tuvsud.com/ps-cert?q=cert.V9.061317.0006.Rev.00</a>				
<b>Report No.:</b>		713210321		
<b>Valid from:</b>		2021-04-01		
<b>Valid until:</b>		2022-05-26		
<b>Date,</b>		2021-04-01		
				
		Christoph Dicks Head of Certification/Notified Body		
Page 1 of 1 TÜV SÜD Product Service GmbH is Notified Body with identification no. 0123				
TÜV SÜD Product Service GmbH • Certification Body • Ridlerstraße 65 • 80339 Munich • Germany				
TUV®				

# Declaration of Conformity

	
<b>Declaration of Conformity</b>	
<b>Manufacturer</b>	Xiamen Boson Biotech Co., Ltd. 90-94 Tianfeng Road, Jimei North Industrial Park, Xiamen, Fujian 361021, P. R. China.
<b>European Representative</b>	Lotus NL B.V. Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.
<b>Product</b>	Rapid SARS-CoV-2 Antigen Test Card, REF: 1N40C5-2 Rapid SARS-CoV-2 Antigen Test Card, REF: 1N40C5-4 Rapid SARS-CoV-2 Antigen Test Card, REF: 1N40C5-6
<b>Classification</b>	Self-Test IVD
<b>Confirmatory Assessment Route</b>	Directive 98/79/EC on In Vitro Diagnostic Medical Devices (IVDD), Annex III section 6
We herewith declare that above mentioned products meet the provisions of the council Directive 98/79/EC for medical devices. All Supporting Documentation is Retained under the premises of the manufacturer. We have sole responsibility for issuing the Declaration of Conformity.	
<b>Standard Applied</b>	ISO 13485:2016                      EN ISO 14971:2012 EN ISO 18113-1:2011            EN ISO 18113-4:2011 ISO 15223-1:2016                EN 62366-1:2015 EN 13612:2002                    EN ISO 23640:2015
<b>Notified Body</b>	TÜV SÜD Product Service GmbH, Ridlerstraße 65, 80339 Munich, Germany. Identification number: 0123
<b>(EC) Certificate(s)</b>	No. V9 061317 0006 Rev. 00
<b>Start of CE-Marking</b>	2021-04-01
<b>Place, Date of Issue</b>	Xiamen, 2021-04-01
<b>Signature</b>	 (Signed By Boson Representative) Name: Changgong Zhang Title: General Manager
	 0123

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# Company Information

## Boson Introduction

Xiamen Boson Biotech Co., Ltd., as a specialist of in vitro diagnostic kits field was founded in 2001, we develop and manufacture high-quality point of care and other immunoassay kits for world-wide market.

Our 10,000 square meter facility is operated strictly under ISO 13485:2016 and GMP guidelines.

Our product lines provide the immunoassays in various formats to detect cardiac markers, drugs of abuse, fertility hormones, infectious diseases, tumor markers and animal diseases. Many of our products have been approved by NMPA and CE.

Our company presents its product with well-designed Boson and HomeScan packaging. We also provide OEM and private label service for customers.

Xiamen Boson Biotech's mission is to provide the affordable high-quality products to help fight diseases and illicit substance abuse.

# Medical Device Production License

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<h3>Medical Apparatus and Instrument Production License</h3>	
License Number: Production License No. 20100174 from the Fujian Drug Administration	
<b>Business Name:</b> Xiamen Boson Biotech Co., Ltd.	<b>Manufacture Address:</b> 90-94 Tianfeng Road, Jimei North Industrial Park, Xiamen
<b>Legal Representative:</b> Changgong Zhang	<b>Manufacture Scope:</b> Class II/III 6840 in-vitro diagnostic reagent 22-04 immunoassay equipment from Revised China Medical Device Classification
<b>Person in Charge:</b> Changgong Zhang	
<b>Address:</b> 90-94 Tianfeng Road, Jimei North Industrial Park, Xiamen	<b>Issuing Authority:</b> Fujian Drug Administration
<b>Expiration Date:</b> November 5, 2024	<b>Issuing Date:</b> December 31, 2019

Prepared by China Food and Drug Administration

## Our Market

Our products have been exported to more than 70 countries: Britain, Germany, Italy, Brazil, Argentina, Australia, India, Malaysia, South Africa, etc.

*Please contact us or your agent for details.*